

# Gynecological Products

IndoSurgicals manufacture and supply superior quality obstetrics and Gynecology Products to satisfy the needs of leading obstetrics and gynecology medical professionals around the world. IndoSurgicals focused on better health outcomes for women, which is why we partner with gynecologist and obstetrics to continually deliver innovative products and the best medical device technology. We offer a full range of high-quality obstetrics and gynecology products including disposables, instruments and equipment.



## Products Inside

- Vacuum Extractor Sets - Manual
- Vacuum Extractor Sets - Electric
- Umbilical Cord Clamp
- Gynecological Aspiration Kit
- Latex Gynecological Gloves
- Breast Pumps
- Vaginal Speculum
- Pessaries



## 55000 - Vacuum Extractor Set, Manual Operated, Malmstorm Type

### Specification:

- Vacuum Extractor Set, Malmstorm type, manual operated.
- Complete with 3 S.S. Cups, (40,50, 60mm)
- PVC tubes
- Traction handle
- Suction pump
- Glass jar in between top and bottom metal plates fitted with vacuum gauge and control valve.
- Packed in case.



## 55001 - Vacuum Extractor Set, Manual Operated

### Specification:

- Complete with 3 Silicon Cups (without valves), (40, 50, 60mm)
- Autoclavable Silicone pressure tube
- Suction pump
- Glass jar fitted between M.S. Basket with vacuum gauge and control valve
- Packed in case



## 55002 - Vacuum Extractor Set, Manual Operated

### Specification:

- Vacuum Extractor Set, Manual operated complete with 3 stainless steel cups (40,50,60mm)
- Autoclavable silicone pressure tube
- Traction handle
- Suction pump
- Glass jar fitted between M.S. Basket with vacuum gauge and control valve
- Packed in case





## 55003 - Vacuum Extractor Set, Manual Operated

### Specification:

- Complete with 3 stainless steel cups (40, 50, 60mm) & 2 Silicon Cups (without valves)
- Autoclavable silicone pressure tube
- Traction handle
- Suction pump
- Glass jar fitted between M.S. Basket with vacuum gauge and control valve
- Packed in case.



## Silicon Cups for Vacuum Extractor

| Item Code | Opening Dia |
|-----------|-------------|
| 55006A    | 50MM        |
| 55006B    | 60MM        |
| 55006C    | 70MM        |



## 55011 - Abortion Set - Manual

### Specification:

- Abortion Set, Manual operated for medical termination of pregnancy
- Complete with 4 cannulae
- Rubber pressure tubes
- Glass jar fitted with gauge and control valve
- Packed in plastic case.





### 55004 - Vacuum Extractor Electric With S.S. Cups

#### Specification:

- Housing: ABS Plastic Moulded
- Capacity: -575mmHg  $\pm$  10 , 45 Ltrs/min
- Pump Type: Diaphragm Pump
- Jars: 1 x 2 Ltrs Polycarbonate Jar
- Filter: Bacterial Filter Autoclavable/Reusable
- Tubing: 8mmID x 2mtr. (Silicone)
- Vacuum Gauge: 6.25 cm dia 0-760mmHg calibration.
- Power: 220 V AC, 50 / 60 Hz, 90Watt
- Noise Level: 55 dB A  $\pm$  3
- Overflow Safety: Mechanical / Electronic Acoustic Alarm
- Supplied with 3 S.S. cups (40, 50, & 60mmdia.)



### 55005 - Delivery Aid for Obstetricians : Safe for Mother/Baby

#### Specification:

- House : Engineering moulded ABS plastic cabinet with built-in handle and jar hanger.
- Pump Type : Oil-free piston pump.
- Capacity : 0-575mm Hg  $\pm$ 10. Fine regulation valve for necessary adjustment with break mechanism under thumb in ventouse cup. Hi/Lo switch provide full or half vacuum with one flick.
- Vacuum Gauge : Bourden type 5 cms. 0-760mm Hg calibration.
- Filter : Bacterial filter
- Power : 220 V 50/60 Hz, Single phase
- Venouse cup : Silicone - 60mm and 70mm : minimizes risk of trauma and cephalohematoma of fetal cranium.



### 55009 - Umbilical Cord Clamp

#### Specification:

- Umbilical Cord Clamp
- Sterile
- Pack of 100 Pcs.
- Openable / Non-Openable





## 55008 - Gynecological Aspiration Kit (MVA Kit)

### Specification:

- Consisting of Syringe Double Valve (Sz. 60cc),
- Set of 8 pcs sterile cannula (size from 4 to 12mm) with Adapter
- Lubricating oil



## 55010 - Latex Gynecological Gloves

### Specification:

- Latex Gynecological, Long cuff Gloves
- Sterile



| Item Code | Size   |
|-----------|--------|
| 55010-S   | Small  |
| 55010-M   | Medium |
| 55010-L   | Large  |

## 55012 - Manual Breast Pump with Silicon Breast Shield

### List of Components:

- 1. Lever Handle 2. Connector 3. Silicon Adaptor shield 4. Funnel 5. 2 - Component Connector 6. Value Head 7. Value Membrane 8. Feeding Bottle 9. Bottle Cover 10. Nipple 11. Bottle Connector 12. Bottle Stand

### Features:

- Pull Handle Design. Its dual-ring silica gel unique 4 air inlet design makes suction around mamma uniformly and alleviate pain.

### Technique Parameter:

- Stimulation phase : 0 ---- 150mmHg  
0 ---- 20kPa
- Expression phase : 0 ---- 150mmHg  
0 ---- 20kPa





## 55013 - Vaginal Speculum - Cusco (Stainless Steel)

| Item Code | Size   |
|-----------|--------|
| 55013S    | Small  |
| 55013M    | Medium |
| 55013L    | Large  |



## 55014 - Vaginal Speculum - Duck Bill (Sims) (Stainless Steel)

| Item Code | Size   |
|-----------|--------|
| 55014S    | Small  |
| 55014M    | Medium |
| 55014L    | Large  |



## 55015 - Vaginal Speculum - Graves (Stainless Steel)

| Item Code | Size   |
|-----------|--------|
| 55015S    | Small  |
| 55015M    | Medium |
| 55015L    | Large  |





Ring Vaginal Pessary Silicone Non Sterile

Specification:

- IndoSurgicals Ring Vaginal Pessary are 100% silicone pessaries that are used in the management of vaginal prolapse of the urogenital organs.
- These pessaries are soft and pliable, making them easy to insert and remove. Available in 3 sizes and all are without support.
- Ring Pessary is very common for a first to second degree prolapse. The Ring with support can also be used on an accompanying cystocele. Insertion is eased with the folding action of the Ring.



| Item Code | Size   | Outer Diameter                          |
|-----------|--------|---|
| 55020     | Small  | 70 mm and Inner Diameter - 50 mm approx |
| 55021     | Medium | 80 mm and Inner Diameter - 60 mm approx |
| 55022     | Large  | 87 mm and Inner Diameter - 67 mm approx |

LABORATORY REAGENTS  
& FINE CHEMICALS



- Dehydrated Culture Media
- Laboratory Reagents
- Analytical Reagents
- Fine Chemicals
- Rare Earth Oxides

- HPLC Solvents
- High Purity Acids
- PH Indicators & Solutions
- Biological Stains
- Concentrated Volumetric Solutions