

HEIGHT & WEIGHT SCALES



Height Measuring Scales

- » Height Measuring Scale, floor model
- » Measuring in inches and centimeters
- » Height stand disassembles into 4 sections & easy to transport in disassembled condition
- » Suitable for height measuring at any location
- » Accurate, portable, durable and light weight with all plastic parts
- » Suitable for medical practices, hospitals, nursing home, school, home etc.
- » Supplied with knock down condition with how to assemble manual in English language

SKU	Type	Measuring range
20015	Without Weighing Scale	20 - 210cm with 1mm Graduation
20016	With Mechanical Weighing Scale 130Kg	20 - 205cm with 1mm Graduation
20017	With Digital Weighing Scale 150Kg	20 - 205cm with 1mm Graduation

Infantometer Baby Measuring Device

- » Sleek broad acrylic base with one sliding side adjustable as per length of baby, dual scale for direct reading in cms / inches upto 90cm
- » Folding sides for easy storage, baby can comfortably lay on the chest of the scale.

SKU	Type
20014	18" x 5" (Length x Width)
20024	18" x 7" (Length x Width)

20019 - Infantometer

- » The IndoSurgicals Infantometer is lightweight and the convenient fold-up mechanism makes it compact and easy to transport.
- » The removable foot stop is securely guided along rails and the measuring boards smooth surface makes it easy to clean.
- » Simple operation for quick mobile use
- » Folds up for space-saving storage High-quality materials for long service life
- » Smooth surfaces for easy cleaning
- » Measuring range 4 to 39 inch

20040 - Infantometer Rod Type

- » Measuring in inches and centimeters
- » Infantometer disassembles into 2 sections & easy to transport in disassembled condition
- » Suitable for height measuring at any location
- » Accurate, portable, durable and light weight with all plastic parts
- » Suitable for medical practices, hospitals, nursing home, school, home etc.



HEIGHT & WEIGHT SCALES



20018 - Stature Meter, Height Measure Scale (Roll-up Model) 200cm

- » Roll-up model for space-saving
- » Display window to read-out of results easily
- » Durable metal tape
- » Hassle free fixing on wall with only one screw, Package Contents Stature Meter and all fixing accessories
- » Measuring range : 0-200 cm, Division value 0.1cm

20030 - Head Circumference Measuring Tape

- » This head circumference measuring tape is with unique 3-slot principle for accurate measurement.
- » Measuring range: Upto 59cm.
- » Material: Flexible, Washable, Non tearable & Nonstretchable PVC.
- » Graduation: 1 mm.

Age	Head Circumference	Age	Head Circumference
At birth	34-35cm	6 Months	42-43cm
2 Months	38cm	1 Year	45-46cm
3 Months	40cm	2 Years	47-48cm
4 Months	41cm	5 Years	50-51cm

20021 - Mid Upper Arm Circumference (MUAC) Tape

- » Mid Upper Arm Circumference Tape - Child.
- » MUAC tapes are predominately used to measure the upper arm circumference of children, helping identify malnutrition.
- » Flexible & Washable Material. Measuring in cm. Graduation - 1mm.
- » Measuring Up To 30+ cm.
- » Made from soft non tearable vinyl sheet.
- » Colour Coded : 3-11.5cm (Red), 11.5-12.5cm (Yellow) 12.5 To 30+cm (Green).

MUAC Tape for Child & Adult - Deluxe Quality

- » This measuring tapes are used to measure the mid upper arm circumference (MUAC) of Adults and Children. MUAC tapes are predominately used to measure the upper arm circumference of children, helping identify malnutrition. The tapes are colored at various cut-off points reflecting severity of malnutrition.
- » Material: Flexible, Washable, Non tearable & Nonstretchable PVC.
- » Printed scale and colour can withstand all commercially available cleaning products.

Technical specifications for CHILD MUAC Tape:


- » Measurement range: Upto 26.5cm.
- » Graduation: 1mm.
- » Color coded: Red (sever) 0-11.5cm; Yellow (moderate) 11.5-12.5cm; green (normal) 13.5-26.5cm

Technical specifications for ADULT MUAC Tape:

- » Measurement range: Upto 45.5cm.
- » Graduation: 1mm.
- » Color coded: Red (sever) 0-21cm; Yellow (moderate) 21-23cm; green (normal) 23-45.5cm.



Double Function : Head Circumference with front side of the tape and facial symmetry the rear side of the tape.



SKU	Type
20035	Child
20036	Adult

HEIGHT & WEIGHT SCALES



EQBR9201 – Equinox Personal Weighing Scale Mechanical (130 kg)

- » The weighing scale has a capacity of up to 130 kg.
- » It is built with durable construction that accommodates all body types.
- » The weighing scale comes with a zero adjuster, which needs to be adjusted before stepping on the scale.
- » Equinox weighing scale EQ-BR-9201 is a solidly constructed weighing scale with a platform size of 26.7 x 26.8 x 5.0 cm.
- » Easy to use without any hassle of installing or adding a battery.
- » Digits are big and clear.



EQBE9300 – Equinox Personal Weighing Scale Digital (150 Kg)

- » Stylish look – Slim and sleek which seamlessly fits into any decor
- » Display – Large LCD display makes it comfortable for all our users to view results on the device
- » Auto-step-on – Product comes with auto-step-on feature which means all you need to do is step on the scale and your readings will reflect on the LCD display
- » Power – CR2032, 3 volts lithium battery
- » Low battery/error/overload indicator
- » Auto off function – Once you step down the device will automatically turns off within a few seconds. The reading is displayed on the LCD display for 5 seconds before it is turned off



20007 – Baby Weighing Scales (Pan Type)

- » Baby Weighing Scales. Plastic Die Cast Body with Plastic Pan.
- » Capacity : 20 Kg. & Graduation : 50 gm.
- » Measuring Tape is provided in the Plastic Pan to measure height of baby.
- » Accurate & Reliable weight. Wide View dial with moving pointer.
- » Easily removable pan.



20008 – Digital Baby Weighing Scale

- » Capacity : 20 Kg. x 5g
- » Height scale to measure length of baby
- » Large scale size : 563 x 333 mm
- » Accurately measures weight in 5g increments up to 20 Kg. With HOLD function for precise weighing of active baby
- » TARE feature subtracts added weight of diaper, clothing or blanket etc.
- » Large display of easy reading, Safety weighing platform
- » Includes 2 x AAA batteries
- » Size : 563 x 333 x 88 mm



20020 – Digital Baby/Child Weighing Machine, Capacity (30 Kg)

- » Range of measurement: 20g-30kg
- » Division: 5g
- » Battery: 4 x 1.5V AA
- » Weight unit: kg
- » Overloaded indication: O_Ld
- » Low power indication: LO
- » Auto switch off, TARE function
- » Working conditions: +0°C + 40°C : RH85%
- » Storage conditions: -10°C + 60°C : RH85%



HEIGHT & WEIGHT SCALES



20100 – Digital Baby Weighing Scales

- » Large LCD Display
- » Max. Capacity : 20 Kg., Graduation : 5 gm.
- » Tare weight function : Hold function for recording weight.
- » Automatic switch off, Battery included



EQBE55 – Equinox Digital Baby Weighing Scale (25 Kg)

- » Large 1.2 inch LCD display
- » Monitors weight and height
- » Built-in memory for 8 infants with 15 records for each infant
- » Portable device
- » Includes delta weight function meaning it will display the difference between the previously recorded weight / height and current weight / height
- » Low battery and overload indicator



Baby Weighing Scale, Salter Type (Dial)

- » With two metal suspension hooks with adjustment screw
- » Capacity: 25 Kg. x 100 Grams

SKU	Type
20010	ABS Plastic Body
20011	Metal Body



20012 – Weighing Trouser for Baby Weighing scales

- » Weighing Trouser for Baby Weighing Scales
- » With strap for suspension from Infant Scale.



Baby Weighing Balance / Scale (Spring type)

- » Capacity: 5 Kg x 100 Gms.
- » Colour Coded with zero adjustment device.
- » Compact to fit in midwifery kits.
- » Light weight – Nylon.
- » With hook, suspension ring & rustproof firm grip handle.

SKU	Capacity
20013	5 kg
20040	10 kg

