

MEDICAL RUBBER PRODUCTS



Air Cushion (Invalid Air Rings)

General: 35000 to 35003 & 35043
Size: 30cm, 35cm, 40cm, 42.5cm & 45cm
Deluxe: 35004 to 35006
Size: 35cm, 40cm & 45cm




Breast Pump

35007: Complete (glass barrel + rubber bulb)
35008: Rubber bulbs for breast pump, 55ml
35009: Rubber bulbs for breast pump, 25ml




Enema Syringe

SKU: 35010
Enema Syringe Higginson, in Red Colour with Rectal Nozzle & Valve



Drainage Sheet Corrugated

SKU: 35011
Size: 150mm X 300mm



Chloroform Bellow

SKU: 35012 & 35013
Size: 500ml & 1000ml



Stomach Pump Tube

SKU: 35024
150 cm with bulb, funnel and plastic mouth gag



Eye/Ear & Ulcer Syringe

SKU: 35014 to 35017
Size: 25ml, 50ml, 75ml & 100ml



Infant / Rectal Syringe

SKU: 35018 to 35020
Size: 25ml, 64ml & 114ml



Vaginal Douche Spray

SKU: 35021
Complete With Plastic Nozzle



Nasal Aspirator

SKU: 35022
Rubber Bulb fitted with Plastic Nozzle



Silicone Nasal Sucker Suction Unit

SKU: 35060
Capacity Approx. 80 ML



Pessary Rubber Ring


SKU: 35023, 35044 & 35045
Size: 50 mm, 75 mm & 87 mm

MEDICAL RUBBER PRODUCTS



Pipette Bulb

Item Code: 35031 to 35033
Size: Small 10ml, Medium 15ml & Large 20ml



Tourniquet (Latex)

Item Code: 35034
Features: Tourniquet Latex Elastic Strip Tubular
Reusable, Non-sterile. Length: 75-100cm Approx.



Tourniquet

Item Code: 35026 & 35027
Color: Blue & Red



Mackintosh Sheet

35028: Deluxe Quality
35029: Regular Quality

- » IndoSurgicals Mackintosh Sheeting (Hospital Sheeting) are made from best quality Indian soft rubber.
- » The sheeting has high breaking and tearing strength.
- » Washability : With cold water and mild soap.
- » Roll length : 10 meter
- » Width approx. : 90cm
- » Colors : Red/Green, Red/Blue.



Douche Bag

Item Code: 35030, 35040 & 35041
Size: 1 Litre, 2 Litre & 4 Litre

- » High Quality - Superior Rubber Latex Bag with Superior Tubing and Smooth Nozzles.
- » Fountain top for easy cleaning & hygiene - Non-leeching & no disturbing odors.
- » No cheap "bottle converted-to-enema" combination systems having leaky adapters - use just for Enemas.
- » Suited for High (Deep) Enemas; Supplied with complete accessories.
- » High quality accessories - hose, nozzles, clamp and a flexible colon tip for deep cleanse.
- » INSTRUCTION MANUAL : Supplied with Instruction Manual in English language.



Hot Water Bottles

Item Code: 35036 to 35038
Color: Red, Blue & Green

- » Made from natural rubber, the flexible bottle conforms to your body
- » Two Side Ribbed
- » There's nothing like an old fashioned hot water bottle to keep you warm at night.
- » Soothe aches and pains by filling the bottle with warm or cold water and placing on the affected area for 10 to 15 minutes per hour.



Kelly's Pad with Pump

Item Code: 35100

- » Kellys Pad is made up of a good quality rubber.
- » Supplied with Pump
- » Kellys Pad is use while delivery.