


# X-Ray Supplies

IndoSurgicals is a manufacturer and global distributor of X-Ray Equipment. Our products are available in more than 50 countries round the globe. Whether you are a reseller, small independent diagnostic center, Radiology, chiropractic or veterinary clinic, or a hospital needing new equipment systems, we are ready to meet your needs.



## Products Inside

- X-Ray View Box
- X-Ray Cassettes
- X-Ray Hangers
- Lead Aprons
- Thyroid Collars
- Lead Gloves
- Lead Goggles
- Gonad & Ovarian Shields
- Lead Apron Racks



## X-Ray View Box, LED, Deluxe Quality with dimmer to adjust brightness of screen



Item Code: 91000



Item Code: 91001

### Specification:

- » Adopt LED light source with life of 50,000 hours
- » Without the effect of temperature difference as CCFL, the viewer can lighted as soon turned on.
- » With high brightness of more than 4000lux. Exclusive for viewing details image freely.
- » Light source, layout in equal spaced points, which is more uniform and the uniformity of luminance is above 90%.
- » Maximum brightness can reach up to 5000cd/m2 which suitable for different density medical films.
- » Individual function for each bank.
- » The light frequency is above 50kHz, effectively relieves the fatigue of vision.
- » The wide power design which fit for the power of different areas, the viewers work well within AC200-250V
- » Ultra thin design and the power set inside the viewers, thickness is aprox. 4.5cm, available for desk and wall

Item Code	91000	91001	91002	91003
	Single Film	Double Film	Triple Film	Four Film
Viewing Area (Inch)	14 x 17	28 x 17	42 x 17	56 x 17
External Dimension (mm)	500 x 512 x 47	850 x 512 x 47	1200 x 512 x 47	1570 x 512 x 47
Color Temperature	>8000K	>8000K	>8000K	>8000K
Light Source	LED	LED	LED	LED
Lamp Life Approx.	100,000 Hours	100,000 Hours	100,000 Hours	100,000 Hours
Power	AC100-220V, 50HZ	AC100-220V, 50HZ	AC100-220V, 50HZ	AC100-220V, 50HZ
Power Supply	45W	90W	135W	180W



## X-Ray View Box, LED, General Quality with dimmer to adjust brightness of screen

### Specification:

- » **Body:** Powder Coated MS Cabinet with Metal Frame
- » **View:** Single, Double, Triple and Four View Available
- » **Sheet:** Acrylic Sheet
- » **Power:** 12 V, 50-60Hz
- » **Bulb:** High Glossy Led Bulbs
- » **Parts:** Heavy Duty Switch & Socket Used
- » **Safety:** Shock Proof Mechanism
- » **Holding:** Plastic Clips for X-ray film Holding
- » **Handling:** Light Weight & Portable



Item Code	Type
91080	X-Ray View Box, LED, General Quality, Single Film
91081	X-Ray View Box, LED, General Quality, Double Film
91082	X-Ray View Box, LED, General Quality, Triple Film
91083	X-Ray View Box, LED, General Quality, Four Film

## X-Ray View Box, Flourescent Tubes

### Specification:

- » Flourescent Tubes, Plastic Body
- » With high brightness



Item Code	Type
91004	X-Ray Viewer Single Film
91005	X-Ray Viewer Double Film
91006	X-Ray Viewer Triple Film
91007	X-Ray Viewer Four Film





## X-ray Cassettes

- Made of specified aluminum alloy frame.
- Specified engineering plastics and high-stamina aluminum plate.
- Without Film / Screen.

Item Code	Type
91012a	6.5" X 8.5"
91012b	8" X 10"
91012c	10" X 12"
91012d	14" X 17"
91012e	12" X 15"
91012f	14" X 14"



## X-Ray Hangers

- Made of Stainless Steel.
- Made from high quality materials.

Item Code	Type
91013a	6.5" X 8.5"
91013b	8" X 10"
91013c	10" X 12"
91013d	14" X 17"
91013e	12" X 15"
91013f	14" X 14"

## Lead Apron (Wrap Around)

### Specification:

- Lead Apron is specially designed for complete protection .
- It gives protection to the whole body.
- Provides maximum comfort, very less fatigue and freedom of movement.
- Front Lead Equivalence : 0.5mm & Back Lead Equivalence : 0.25mm

Item Code	Lead Equivalency
91008	0.35mm
91108	0.50mm





## Lead Apron (Velcro Type)

### Specification:

- It is designed to provide efficacious shield against radiation during X-ray process.
- Provides maximum comfort, very less fatigue and freedom of movement
- With Pocket, Unisex size
- Supplied with Hanger

Item Code	Lead Equivalency
91009	0.35mm
91109	0.50mm



## 91110 - X Ray Lead Apron, Strap Type, Lead Equivalency 0.35mm

### Specification:

- It is designed to provide efficacious shield against radiation during X-ray process.
- Provides maximum comfort, very less fatigue and freedom of movement
- With Pocket
- Unisex size
- Supplied with Hanger



## 91073 - Half Lead Apron

### Specification:

- Half Lead Apron
- Front Lead Equivalence : 0.5mm
- Size Approx. : 55 x 45cm with adjustable strap with buckle
- Weight: Approx. 2 Kg.
- Color: Blue





## Thyroid Collars

### Specification:

- Provide protection for the neck and thyroid area.
- It helps to reduce patient exposure during intraoral radiography.

Item Code	Type
91010	Thyroid Collar
91100	Thyroid Collar, Imported



## 91011 - Lead Gloves

### Specification:

- Radiation resistant gloves
- Provides maximum comfort



## X-Ray Lead Goggles

### Specification:

- Specially designed to fit over the eyes, leaded fitovers are an excellent choice for those requiring good splash protection from radiation.

Item Code	Type
91050	X-Ray Lead Goggles, Wrap Around
91051	X-Ray Lead Goggles





### 91071 - X-Ray Ovarian Shield

#### Specification:

- Provides protection to the Ovaries during procedures where the patient needs ovary protection.
- Lead equivalency is The 0.5 mm.



### 91072 - X-Ray Protection Gonad Shield

#### Specification:

- Triangular style with ties.
- Lead equivalency is The 0.5 mm.
- Specifically to protect the genital area.
- The gonad shield comes with straps to hold it in place.



### 91070 - Lead Apron Trolley

#### Specification:

- Stand made of stainless steel
- Mobile on four wheel castors
- Supplied with five movable hangers for Lead Apron
- Easy to clean

